Pancreatic Cancer Organoids from the Human Cancer Models Initiative Biobank Reflect Disease Genotypes, Capture Patient Heterogeneity, and are Amenable to Therapeutic Screening

Stephen Friend, MS; Matthew Graziano, BS; Rachel Demerath, MS; and James Clinton, PhD ATCC[®], Manassas, VA 20110

Introduction

- Pancreatic cancer is the fourth leading cause of cancerrelated death in the United States (~3% of all cancer cases) and has one of the lowest 5-year relative survival rates.
- There is a lack of clinically representative, easily available, validated, in vitro pancreatic cancer models that reflect the genomic and phenotypic diversity of the disease.
- The Human Cancer Models Initiative (HCMI) is an collaborative project dev international devoted to development and distribution of approximately 1,000 novel human primary tissue-derived tumor models supported with clinical and molecular annotation.
- ATCC[®] has made 300+ of these next-generation cancer models, including 50+ pancreatic cancer organoids (PCOs), available to the research community through our catalog.
- Here, we describe a subset of the pancreatic model cohort, highlight their clinical characteristics, characterize their genotype, and investigate their response to a 10-compound panel of chemotherapeutics that include taxanes, platinum compounds, KRAS-targeting inhibitors, and PARP inhibitors.

Human Cancer Models Initiative (HCMI)

HCMI is an international consortium of organizations with the shared goal of creating ne generation in vitro cancer models that better represent the diversity and complexity of hun cancers than seen in existing cancer cell lines. These models are annotated with detail clinical (e.g., patient demographics and treatment history) and molecular data (e.g., WC WES). These novel models are manufactured and distributed by ATCC[®]. Over 300 mod are currently available, including 3-D cancer organoids, from 23 different primary tissue sit Over 50 PCOs are currently available, a subset of which are described in Table 1.

Methods

Organoid culture: PCOs were sourced from ATCC[®] and cultured in standard extracellular matrix (ECM)-embedded conditions according to recommended protocols (Clinton et al. 2019) utilizing Organoid Growth Kit 1B (ATCC[®] ACS-7101[™]), Cell Basement Membrane (ATCC[®] ACS-3035[™]), and ROCK Inhibitor Y27632 (ATCC[®] ACS-3030™)

Drug response: Viability studies were performed between passages 2 and 10. Passage occurred 72 hours prior to seeding to allow single cells and cell fragments to reform small organoids prior to seeding and dosing. Organoids were collected, ECM removed, and seeded

ATCC



Figure 1: Summary of workflow. Preparation dosing, and assay timepoints.

as intact organoids into white walled 96-well plates at the equivalent of ~2.5x10³ cells/well in complete growth media supplemented with 0.5 mg/mL ECM. Compounds for the drug screening panel were KRpep-2d, MRTX1133, nab-paclitaxel, cisplatin, 5-fluorouracil (5-FU), olaparib, oxaliplatin, PJ34, paclitaxel, and rucaparib. Drugs were reconstituted in either H₂O, DMSO, or NaCI. Positive control for drug response was 200 µM cisplatin. Drug exposure and assay timepoints are shown in Figure 1. Viability was determined using the CellTiter-Glo[®] 3D Luminescent Cell Viability Assay (Promega[®]). Response was normalized to the vehicle condition and expressed as percent viability. Figures were plotted, non-linear curves generated, and IC₅₀ was calculated in GraphPad Prism[®]. Points reflect average with SEM of 5-6 technical replicates. Z' Factor scores were determined using vehicle control and positive control.

Imaging and histology: Organoids in culture were routinely monitored via brightfield microscopy. Prior to each viability assay, images were taken of each well to monitor morphology. For histology, organoids were collected, ECM removed, pelleted, fixed in 4% paraformaldehyde, embedded in paraffin, sectioned at 10 µm, mounted on slides, and stained with hematoxylin and eosin (H&E) to visualize morphology (Figure 3).

Sequencing: Organoids were collected, ECM removed, and DNA extracted using a QIAGEN® EZ1® Advanced XL instrument with an EZ1[®] DNA Tissue Kit. DNA concentration and purity were measured using a NanoDrop[®] 2000. Libraries were prepared using a targeted cancer panel (~600 genes) for hybrid capture. Sequencing was performed by Illumina[®] NovaSeq[®] 6000. Average target depth was >1300X. Low quality data were filtered, and alignment was performed against GRCh37. Variants were filtered and annotated. Oncoplot of selected key genes and variants associated with pancreatic cancer are shown in Figure 2.



10801 University Boulevard, Manassas, Virginia 20110-2209

© 2024 American Type Culture Collection. The ATCC trademarks owned by the American Type Culture Collection unless indicated otherwise. CellTiter-Glo and Promega are registered trademarks of Promega Corporation. GraphPad Prism is a registered trademark of GraphPad Software, LLC. QIAGEN and EZ1 are registered trademarks of the QIAGEN and EZ1 are registered trademarks of Illumina, Inc.

Phone: 800.638.6597

Results

Pancreatic

ancreatic

ancreatio

ancreatic

ancreatic

Pancreatic

ancreatio

Adeno. ductal type

Adeno, ductal type

Adeno. ductal type

Adeno. ductal type

Adeno, ductal type

Adeno. ductal type

deno. ductal type

Adeno. ductal type

Adeno. (NOS)

Adeno. (NOS)

PDM-38[™] DM-289™

PDM-101™

PDM-110™

PDM-134™

PDM-198™

PDM-200[™]

PDM-203™

PDM-204™

PDM-205™

PDM-24™

PDM-270[™]

PDM-28™

PDM-29™

Histological subtype

Anatomic location

Vital status

4% 17%

4% 13%

7%

7%

20%

TCGA HCM

the

	TP53	5 <mark>3%</mark>	74%
	SMAD4	19 <mark>%</mark>	22 <mark>%</mark>
	CDKN2A	15%	65%
	RNF43	6%	4%
ext-	TGFBR2	4%	17%
nan	ARID1A	4%	13%
iled	GNAS	4%	7%
GS,	KDM6A	4%	7%
dels	ATM	3%	15% <mark>/</mark>
tes.	APC	2%	20 <mark>%</mark>

BRCA2



168h



Figure 3: Heterogeneity of morphology across PCOs (n=4). Brightfield and H&E staining reveal typical morphology of PCOs during routine culture, exhibiting either a cystic appearance with defined lumens or without a central lumen. Last column shows impact of toxicity from positive control (200 µM cisplatin) as shown by a loss of structure and cellular organization.

Log [Nab-paclitaxel]

Table 1: Clinical characteristics of established human PCOs (n=46) from the HCMI biobank. Models were predominantly derived from primary tumors (75%), had an adenocarcinoma ductal histological subtype (PDAC, 74%), and were from patients with an average age of 65 years at diagnosis that were primarily female (56%). Highlighted models (n=4) were used in subsequent experiments.

® No.	Cancer Type	Histological Subtype	Туре	Acquisition Site	Gender	Age	Clinical Stage
0 ™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Male	75	Stage III
1™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Female	60	Stage III
2™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Male	57	Stage III
3™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Male	77	
4™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Female	73	
5™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Male	75	
9™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Male	69	Stage III
0™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Female	71	Stage III
1™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Female	61	Stage III
21™	Pancreatic	Adeno. ductal type	Primary	Pancreatic head	Male	54	Stage III
23™	Pancreatic	Adeno. ductal type	Primary	Pancreatic body	Female	67	Stage III
0 ™	Pancreatic	Adeno. (NOS)	Primary	Pancreas	Female		Stage IIA
06™	Pancreatic	Other	Metastasis	Liver	Female	63	Stage IIB
07™	Pancreatic	Adeno. ductal type	Metastasis	Liver	Female	49	Stage IIB
08™	Pancreatic	Other	Metastasis	Liver	Female	53	
26™	Pancreatic	Adeno. (NOS)	Metastasis	Pleural cavity	Male	66	Stage IV
64™	Pancreatic	Adeno. ductal type	Metastasis	Peritoneum	Male	66	
68™	Pancreatic	Adeno. ductal type	Metastasis	Other	Female	60	
70™	Pancreatic	Adeno. ductal type	Metastasis	Pleural cavity	Male	73	
79™	Pancreatic	Adeno. ductal type	Metastasis	Liver	Female	74	
21™	Pancreatic	Infiltrating ductal carcinoma	Metastasis	Lymph node(s)	Male	72	
22™	Pancreatic	Infiltrating ductal carcinoma	Metastasis	Pancreatic head	Male	62	
88 TM	Pancreatic	Moderately differentiated Adena	Metastasis	Liver	Female	68	Stage IV



Histological subtype

ocarcinoma (NC Adenocarcinoma (ductal) filtrating ductal carcinon gnet ring cell carcinom d. differentiated adeno





Run 1	(
Run 2	(
Run 3	(
Run 4	(
Run 5	(
Run 6	(
Average	(

Summary and Conclusions

PDM-289™

Acknowledgements and References

We would like to acknowledge all founders and contributors to the HCMI: The National Cancer Institute, the National Institutes of Health. Cancer Research UK, Wellcome Sanger Institute, foundation Hubrecht Organoid Technology, Broad Institute, Cold Spring Harbor Laboratory, Stanford University, Weill Cornell Medical College, University of Verona, and the Hubrecht Institute. The models and some data are from the HCMI: www.cancer.gov/ccg/research/functional-genomics/hcmi; dbGaP accession number phs001486. • Clinton JM and McWilliams-Koeppen P. Initiation, Expansion, and Cryopreservation of Human Primary Tissue-Derived Normal and Diseased Organoids in Embedded Three-Dimensional Culture. Curr Protoc Cell Bio 82(1): e66, 2019. Hallin J, et al. Anti-tumor efficacy of a potent and selective non-covalent KRASG12D inhibitor. Nat Med 28(10): 2171-2182, 2022.



Email: SalesRep@atcc.org

Web: www.atcc.org

38™ (●), PDM-289™ (●)



bv KRAS

exhibited

The HCMI biobank contains organoids from a diverse collection of patients and disease indications, including from pancreatic cancer.

2.03 5.58 0.82 >10 12.63 >200

PCOs from the HCMI have canonical mutations in key genes (KRAS, TP53, SMAD4, and CDKN2A) seen in patient populations and existing large cancer datasets such as the TCGA.

We used PCOs to validate a viability screening assay with a panel of 10 anti-cancer drugs.

PCOs exhibited variable drug toxicity that may be in part a consequence of genotype.

Learn more at www.atcc.org/hcmi & www.atcc.org/organoids